



Pronex Pneumatic Cervical Traction

For the relief of chronic cervical pain

Help patients successfully manage their chronic condition at home. Pronex is the only traction device that creates an even distraction in the anterior and posterior cervical discs while supporting the cervical curve.

Support and versatility

Pronex cradles the reclining patient's head and neck on two soft foam cushions. One cushion supports the back of the head and the other rests against the top of the shoulders. An air-inflated bellows between them provides up to 35 lbs. of continuously adjustable traction. As the bellows expands, it lifts the head upward and supports the cervical curve, maintaining an even force that creates space between the cervical vertebrae. Lateral neck flexion of 15 and 22 degrees can be obtained utilizing an included flexion wedge. Patients have total control over the amount of traction applied. Squeezing an inflator bulb increases the pressure; a release knob gently reduces it. There is no strain on the jaw's temporomandibular joint (TMJ).

Comfort and simplicity equal compliance

Pronex is portable, comfortable and convenient. There is no assembly, no weights, no cables, no springs or levers. It is also restraint free. All of these features combine to encourage regular use by the patient at home, helping to break the chronic pain cycle.



The Pronex System
includes a lightweight carrying case,
headstrap, flexion wedge, instructional
DVD and Operation Manual.

Pronex Features

• Patented design that provides anatomically correct traction

• User friendly and without restraints

• Three unisex sizes:

<u>Size</u>	<u>Measurement</u>	
Regular	14 to 16"	(Fits most women of typical height, weight and neck size)
Large	16 to 18"	(Fits most men of typical height, weight and neck size)
Wide	18" and up	(Fits individuals with large shoulders and neck, or a short neck)

• Approved for HCPCS Code: E0855

Benefits

- Comfortably supports head and neck
- Exerts force evenly and gently
- Applies no pressure on the temporomandibular joint
- Supports normal curvature of cervical spine
- Gentle traction release
- Provides up to 35 lbs. of continuously adjustable traction
- Compact and portable — encourages a high rate of patient compliance

- No assembly, no weights, no cables, no springs or levers
- Encourages regular use

- Reimbursement from Medicare and private insurers is dependable
- Description: Cervical traction equipment not requiring additional stand or frame

PRESCRIBING INFORMATION

RECOMMENDED FOR USE:

The Pronex Cervical Traction Device is intended to treat a musculoskeletal or neurologic impairment of the cervical spine to relieve pain, relax muscle spasms, and decompress spinal structures.

NOT RECOMMENDED FOR USE:

Cervical Traction should be avoided in any condition of the cervical spine where movement can aggravate the condition or result in spinal instability, spinal injury and/or nerve root injury at risk for causing paralysis or ischemia.

CAUTION: Federal law restricts sale or use of this device on the order of a licensed medical professional, under the laws of the state where licensed to practice. The Pronex Cervical Traction Device should only be used as instructed by a medical professional and it should only be used by the person for whom it was prescribed.
THIS PRODUCT CONTAINS NATURAL RUBBER LATEX WHICH MAY CAUSE ALLERGIC REACTIONS.

Please review the Operations Manual provided with the Pronex Pneumatic Traction.

IMPORTANT: Traction should never cause pain. A stretching sensation or some slight discomfort may be felt during treatment. If pain is increased in the neck, back, arms or legs, or if you experience dizziness, nausea or any other type of pain or discomfort during or after treatment, discontinue use immediately and consult your health care provider before further use.

The Pronex Pneumatic Cervical Traction is covered by one or more of U.S. Patent No.(s) 5,916,185; 5,454,0781; 5,441,479. Additional applications pending.